

ELISA Protocol

Wash Buffer:

	PBS
0.1%	Tween-20

Blocking Buffer:

	PBS
0.1%	Tween-20
3%	BSA

AP Base Buffer:

96ml	di-Ethanolamine
1180 μ l	2M MgCl ₂
800ml	H ₂ O

- 1 Coat antigen (5 μ g/ml in PBS) in 150 μ l in Costar EIA/RIA flat bottom 96 well plate (No. 9017).
- 2 Cover the plate and incubate @ 4C overnight.
- 3 Decant by flicking the plate and incubate 30mins. w/ 150 μ l blocking buffer @ 37C. Wash 3x w/ wash buffer.
- 4 Fill the wells w/ 100 μ l of serum samples diluted in wash buffer (start w/ 1:10), cover and incubate 1h @ RT.
- 5 Flick the plate and wash 3x w/ wash buffer.
- 6 Fill the well w/ 100 μ l anti-mouse IgG-AP (1/1000-1/2500 in wash buffer), cover, and incubate 1h @ RT.
- 7 Flick the plate and wash 3x w/ wash buffer.
- 8 Prepare substrate solution by dissolving one Phosphatase Substrate Tablet (Sigma No. 104-105) in 5ml AP Base Buffer (=> 1mg/ml).
- 9 Fill the wells w/ 100 μ l of substrate solution and read the color development in a microplate reader @405nm. The reaction can be stopped by adding 100 μ l 1M NaOH and preserved covered @ -20C.